

# ZEISS Smart Services Dashboard

Monitor and analyze machine data  
in real-time



Seeing beyond

Overview and comparison of the wide range of functions of ZEISS Smart Services Dashboard PLUS compared to the ZEISS Smart Services Dashboard FREE version.

## ZEISS Smart Services Dashboard FREE vs. PLUS

<b>Main features</b>	<b>ZEISS Smart Services Dashboard FREE</b> Real-time monitoring of each ZEISS device	<b>ZEISS Smart Services Dashboard PLUS</b> Real-time monitoring and analysis of all ZEISS machines devices individually and collectively
<b>Analyses over Time – Utilization</b>	Data evaluation of the last 24 hours or 7 days for one measuring system	Unlimited and customized data evaluation of one or more measuring systems or of a complete machine park for comparison of the devices
<b>Analyses over Time – Collisions</b>	Data evaluation of the last 24 hours or 7 days for one measuring system	Unlimited and customized data evaluation of one or more measuring systems or of a complete machine park for comparison of the devices
<b>Analyzes According to Shifts</b> utilizations, collisions, measurement plans and measurements	Data evaluation of one shift only	Data evaluation of customized shift settings
<b>Data Export</b>	Not possible	Data export as .csv file for: utilization, usage, collisions, measurement plans and temperatures
<b>Machine Usage</b> current measurement plan, utilization, machine usage per day, usage protocol	Data evaluation of the last 24 hours or 7 days	Unlimited and customized data specifications
<b>Measurement Capability</b> Temperature monitoring: workpiece, pinhole and CPU ...	Data evaluation of the last 24 hours or 7 days	Unlimited and customized data specifications
<b>Measurement Capability</b> Collision monitoring: 2D and 3D displays	Data evaluation of the last 24 hours or 7 days	Unlimited and customized data specifications
<b>Operating Data</b>	Data evaluation of the last 24 hours or 7 days	Unlimited and customized data specifications
<b>Machine Logbook</b>	Unlimited availability	Unlimited availability
<b>General Machine Data</b>	Unlimited availability	Unlimited availability
<b>Contacting ZEISS</b>	Via the ZEISS support hotline or creation of a support ticket	Via the ZEISS support hotline or creation of a support ticket
<b>Notifications / Alerts</b>	Via browser and e-mail	Via browser and e-mail